



BestFit™ Nasal Mask

Designed for maximum comfort with an advanced cushion design allowing active control of the seal.

BestFit™ Full Face Mask with Anti-asphyxia Valve

Easy to clean at home and sterilizable for multipatient use in a hospital or institutional environment.



Well polished silicone cushion for optimal comfort and fit

Light headgear with breathable, comfortable and durable fabric

Unique feeding tube access port on full face masks

Anti-asphyxia valve closed at 2cm H₂O for superior breath detection

REFERENCE / DESCRIPTION

100010	BestFit Nasal Mask, vented, L
100020	BestFit Nasal Mask, vented, M
100030	BestFit Nasal Mask, vented, S
100040	BestFit Nasal Mask, vented, MW
101010	BestFit Nasal Mask, non-vented, L
101020	BestFit Nasal Mask, non-vented, M
101030	BestFit Nasal Mask, non-vented, S
101040	BestFit Nasal Mask, non-vented, MW
200010	BestFit Full Face Mask, vented, L
200020	BestFit Full Face Mask, vented, M
200030	BestFit Full Face Mask, vented, S
201010	BestFit Full Face Mask, non-vented, L
201020	BestFit Full Face Mask, non-vented, M
201030	BestFit Full Face Mask, non-vented, S
200110	BestFit Full Face Mask, vented w/ AAV, L
200120	BestFit Full Face Mask, vented w/ AAV, M
200130	BestFit Full Face Mask, vented w/ AAV, S

CE certification CE 0123

